

Medical Spill Kits



Cytotoxic
Bio-Waste
Laboratory
Chemical
Veterinary



Contents

Cytotoxic Body Fluids Spill Kit	3
Cytotoxic Chemical Spill Kit	4
LSL "SINGLE" Bio-Waste Spill Kit	5
LSL "COMPACT" Bio-Waste Spill Kit	6
LSL "RESPONSE BIN" Bio-Waste Spill Kit	7
Liquid Spill Lock - Compact Carry Bag Laboratory Spill Kit	8
Liquid Spill Lock Formaldehyde Spill Kit	9
Liquid Spill Lock - Mercury Spill Kit	10
LSL Compact Sharps Collection Kit	11

SDS for Spill Kits featured in this brochure are available on request.



Cytotoxic Body Fluids Spill Kit

Body Fluids Spill Kit

A well designed kit aimed at general clean up of all body fluids dispersed from a patient receiving Chemotherapy treatment. It greatly reduces the possibility of contamination by providing leading edge technology absorbents, as well as the complete personal protection equipment needed in cleaning the whole area up and ensuring the contaminated surround is made aware to others in the area.

The leading edge LSL absorbent offers many features and benefits for "dry cleaning" up cytotoxic body fluid spills:-

- Natural, non-toxic absorbent powder.
- Rapidly immobilizes and deodorizes potentially infectious spills such as blood, vomit, faeces and urine.
- Safe, easy removal, minimizing the risk of personal contact.
- Helps avoid cross contamination.
- Encapsulates solids.
- Unique absorbent properties solidifies all liquid/chemical spills on contact.
- Captures vapor release.
- Chemosorb pad for absorbing any cytotoxic fluid.

For added protection our kit has no sharp edges - our scoop and scraper has rounded edges to avoid puncturing of protective gloves. Also included is a vomit bag for convenience, red alginate bag for the patients contaminated clothing and all the personal protection equipment needed when cleaning up the spill.



Cytotoxic Body Fluids Spill Kit

Product Code: 25261

Qty	Items
1	Instructions
1	Long sleeve impervious gown
1	Face mask with eye shield
1	Hair net
2	Pairs Nitrile gloves
1	500g LSL absorbent powder
1	Chemosorb pad
1	Vomit bag
2	Purple Cytotoxic waste bags
1	Red alginate bag
2	Cable ties
1	100ml Cavicide disinfectant / cleaner (TGA Approved)
2	Cytotoxic spill signs
1	Disposable scoop & scraper
6	Cleaning cloths

For the effective clean up of:



Blood



Vomit



Faeces



Urine



Cytotoxic Chemical Spill Kit

Chemical (Liquid & Powder) Spill Kit

Cytotoxic and other cancer fighting drugs create a hazard to healthcare workers on a daily basis. With this in mind, we have a specific kit designed for the safe and speedy clean up of Cytotoxic chemical based spills. This kit includes everything needed to ensure safety to the designated person cleaning up the spill including full body personal protection equipment, our world class absorbent that "Rapidly Wicks Up" and draws in the liquid spill and the equipment needed to clean and disinfect the spilt area.

This kit can also be used to clean up body fluids; however, we suggest the use of the designated and cost effective Cytotoxic Body Fluid Spill Kit.

The world class LSL absorbent offers many features and benefits for "dry cleaning" cytotoxic chemical spills:-

- Natural, non-toxic absorbent powder.
- Rapidly immobilizes and deodorizes potentially infectious spills such as blood, vomit, faeces and urine.
- Unique absorbent properties solidifies all liquid/chemical spills on contact.
- Safe, easy removal, minimizing the risk of personal contact.
- Helps avoid cross contamination.
- Captures vapor release.
- Chemosorb pads for absorbing liquid cytotoxic fluids and managing powder spills.

Along the same lines as the Cytotoxic Body Fluid Spill Kit and for added protection we ensure that the kit has no sharp edges, our scoop and scraper have rounded edges to avoid puncturing of protective gloves.



Cytotoxic Chemical Fluids Spill Kit

Product Code: 25261CHEM

Qty	Items
1	Instructions
1	Long sleeve impervious gown
1	P2 flat fold mask
1	Hair net
4	Pairs Cytotoxic gloves
1	500g LSL absorbent powder
4	Chemosorb pads
2	Purple Cytotoxic waste bags
1	Red alginate bag
2	Cable ties
1	Bottle of water
2	Cytotoxic spill signs
1	Disposable scoop & scraper
6	Cleaning cloths
1	Pair impervious over boots
1	Protective eyewear
1	100ml Cavicide disinfectant / cleaner (TGA Approved)

For the effective clean up of:



Blood



Vomit



Faeces



Urine



1.5L SINGLE USE Bio-Waste Spill Kit

General Body Fluids Spill Kit

This Single Use Bio-Waste kit offers many features and benefits as the main absorbent is a Natural Mineral - Non-Toxic absorbent Powder which rapidly immobilizes and deodorizes potentially infectious spills. This kit has unique absorbent properties that solidifies spills on contact. The Bio-Waste kit makes for safe and easy removal of the spill contaminate, leaving the effected area clean and disinfected. The kits are ideal for areas in Hospitals, Nursing Homes, Schools, Public Transport, Taxis, Doctors Surgery, Day Care Centers, Restaurants, Pubs & Clubs, Automobiles, Cruise Ships, Veterinary Clinics and Waste Management.

The leading edge LSL absorbent offers many features and benefits for cleaning up spills:-

- Natural, non-toxic absorbent powder.
- Rapidly immobilizes and deodorizes potentially infectious spills such as blood, vomit, faeces and urine.
- Unique absorbent properties solidifies spills on contact.
- Safe, easy removal, minimizing the risk of personal contact.
- Helps avoid cross contamination.
- Encapsulates solids.

For added protection our kit has no sharp edges - our scoop and scraper have rounded edges to avoid puncturing of protective gloves, and all the personal protection needed when cleaning up the spill. Designed as a single use kit, see our "Big Bucket" and "Response Bin" for larger volume spills.



SINGLE USE Bio Waste Spill Kit

Product Code: 25851.5BW

Qty	Items
1	Pair Disposable Nitrile Powder Free Gloves
1	Disposable Apron, One Size fits all
1	Face Mask & Eye Shield
1	Pair Disposable CPE Elastic Shoe Covers
1	TGA Hospital Grade Disinfectant Spray, 100mL
1	Bio-Waste Infectious Waste Absorbent, Sachet, 100g
1	Cloth Wipes
1	Clinical Waste Contaminated Waste Disposal Bag
1	Disposable Scoop & Scraper
1	Instructions plus SDS for TGA Hospital Grade Disinfectant

Optional Extras

Product Code: 2565.5BW

- LSL Bio-Waste 500g Sachet



Product Code: 2565J5.BW

- LSL Bio-Waste 5.5ltr Jug



For the effective clean up of:



Blood



Vomit



Faeces



Urine

15L BIG BUCKET Bio-Waste Spill Kit

General Body Fluids Spill Kit - Large

The BIG BUCKET Bio-Waste Spill kit offers many features and benefits as the main absorbent is a Natural, Non-Toxic absorbent Powder which rapidly immobilizes and de-odorizes potentially infectious spills. This kit has unique absorbent properties that solidifies spills on contact. The Bio-Waste kit makes for the safe and easy removal of the spill contaminates, leaving the effected area clean and disinfected. The kits are ideal for use in Hospitals, Nursing Homes, Schools, Public Transport, Doctors Surgery, Day Care Centers, Restaurants, Pubs & Clubs, Cruise Ships, Veterinary Clinics and Waste Management.

The leading edge LSL absorbent offers many features and benefits for cleaning up spills:-

- Natural, non-toxic absorbent powder.
- Rapidly immobilizes and deodorizes potentially infectious spills such as blood, vomit, faeces and urine.
- Unique absorbent properties solidifies spills on contact.
- Safe, easy removal, minimizing the risk of personal contact.
- Helps avoid cross contamination.
- Encapsulates solids.

For added protection our kit has no sharp edges, our scoop and scraper have rounded edges to avoid puncturing of protective gloves, and all the personal protection needed when cleaning up the spill. Consider our "Handy" kit for single use, easy to store and carry for small spills and the "Response Bin" for larger spills.



COMPACT Bio-Waste Spill Kit

Product Code: 25820BW

Qty	Items
5	Pair Disposable Nitrile Powder Free Gloves
5	Disposable Apron, One Size fits all
5	Face Mask & Eye Shield
5	Pair Disposable CPE Elastic Shoe Covers
5	TGA Hospital Grade Disinfectant Spray, 100mL
1	Bio-Waste Infectious Waste Absorbent, Sachet, 500g
5	Cloth Wipes
5	Clinical Waste Contaminated Waste Disposal Bag
5	Plastic Disposable Scoop & Spatula
1	Instructions plus SDS for TGA Hospital Grade Disinfectant

Optional Extras

Product Code: 232015

- Stronger Virkon disinfectant



For the effective clean up of:



Blood



Vomit



Faeces



Urine

LSL "RESPONSE BIN" Bio-Waste Spill Kit

General Body Fluids Spill Kit - Extra Large

The "Response Bin" Bio-Waste kit offers many features and benefits as the main absorbent is a Natural, Non-Toxic absorbent Powder which rapidly immobilizes and deodorizes potentially infectious spills. This kit has unique absorbent properties that solidifies spills on contact. The Bio-Waste kit makes for a safe and easy removal of the spill contaminants, leaving the effected area clean and disinfected.

The kits are ideal for use in large spill ares such as Hospitals, Correctional Services, Nursing Homes, Veterinary Clinics, Day Care Centers, Schools, Public Transport, Doctors Surgery, Restaurants, Pubs & Clubs, Cruise Ships and Waste Management.

The leading edge LSL absorbent offers many features and benefits for cleaning up spills:-

- Natural, non-toxic absorbent powder.
- Rapidly immobilizes and deodorizes potentially infectious spills such as blood, vomit, faeces and urine.
- Unique absorbent properties solidifies spills on contact.
- Safe, easy removal, minimizing the risk of personal contact
- Helps avoid cross contamination.
- Encapsulates solids.

The wheelie bin makes it easy to transport the kit from area to area when spills occur.

For the effective clean up of:



Blood



Vomit



Faeces



Urine



**SPILL
CONTROL
STATION**

120L LSL "RESPONSE BIN" Bio-Waste Spill Kit

Product Code: 258120BW	
Qty	Items
1	Instructions
1	30L LSL absorbent powder
1	500ml Cavicide disinfectant spray bottle
2	Pairs protective gloves
1	Pair disposable overalls
1	Pair protective goggles
5	Pairs disposable boot covers
2	Contaminated waste bags & cable ties
1	Security lock tag
1	Squeegee
1	Broom
1	Dustpan & brush
1	Brush & scoop

240L LSL "RESPONSE BIN" Bio-Waste Spill Kit

Product Code: 258240BW	
Qty	Items
1	Instructions
2	30L LSL absorbent powder
1	500ml Cavicide disinfectant spray bottle
2	Pairs protective gloves
2	Pair disposable overalls
2	Pair protective goggles
5	Pairs disposable boot covers
10	Contaminated waste bags & cable ties
1	Security lock tag
1	Squeegee
1	Broom
1	Dustpan & brush
1	Brush & scoop

Liquid Spill Lock - Compact Laboratory Spill Kit

Manage your spills of ALL fluids in the lab with our...

High Performance "Uni-Safe pH Plus" Absorbent Binder for all Chemicals and Fluids

Features:

- "Just one absorbent for all fluids".
- Binds all liquids into a non flowing gel.
- High absorbent capacity.
- Acts as a fume/emission barrier – suppresses vapors.
- Changes colour to indicate pH range.
- No hazardous reactions with chemicals.
- Non Combustible – Non Toxic – Not Dangerous Goods.

UNI-SAFE is a powder that absorbs and gels all chemicals, biotech liquids and even body fluids, without creating a reaction. Other absorbents on the market are chemical-specific, so you need several different products. This creates the risk of using the wrong absorbent on the chemical spill or the absorbent MSDS indicates restrictions on some chemicals, creating a potentially dangerous reaction.

Since UNI-SAFE is universal, you save money by only needing one absorbent in your spill kit. UNI-SAFE is a chemical and oil absorbent binder medium universally applicable to all liquids. UNI-SAFE has been developed for and by the German Chemical Industry, in co-operation with the German Fire Brigades, to combat all hazardous accidents. Its material composition will not cause any hazardous reaction not even with strongly reacting and oxidizing substances, such as nitric acid, hydrofluoric acid, and chloresulfonic spills. It will also absorb bromine, sulphuric acid, hydrogen peroxide and caustics without problems. UNI-SAFE acts as an emission blocker. It creates a 'fume barrier' protecting personnel and environment, even in the case of oleum and fuming chloresulfonic acid will stop immediately. UNI-SAFE binds most liquids into a non flowing gel.


Cleaning up spills made safe & easy.



Liquid Spill Lock Compact Carry Bag Laboratory Spill Kit

Product Code: 258510LAB

Qty	Items
1	Instructions
1	Compact carry bag
1	500ml Uni-Safe PH plus binder flip top shaker bottle
2	Yellow Hazchem absorbent pads
1	Pair protective gloves
1	Half face respirator mask & chemical filter ABEK1
1	Chemical splash goggles
1	Waste bag & cable tie
1	Disposable scoop & scraper



Special Feature....

The above chart supplied on each Uni-Safe pH Plus container, shows the indicator range of pH 1 to pH13 which assists with the identification of the type of chemical spilt or to be controlled. Immediate result once Uni-Safe contacts the fluid.

Liquid Spill Lock - Formaldehyde Spill Kit

A safer way to deal with infectious spills

LIQUID SPILL LOCK FORMALDEHYDE NEUTRALISER - Polyform-FTM is the spill control product of choice for the safe and rapid control of harmful vapors from Formaldehyde, Glutaraldehyde, 10% buffered Formalin, and other aldehydes, including Formaldehyde/Phenol* based embalming solutions.

(The Formaldehyde is neutralised, irreversibly polymerised and phenol vapour is subsequently reduced in this reaction by more than 90%)

Once Polyform-FTM free flowing granules are applied, the harmful vapors are eliminated in 2 to 3 minutes. The formaldehyde and other aldehyde solution are neutralised and irreversibly polymerised. This reaction leaves a non-hazardous biodegradable polymer, safe and simple for personnel to clean up and dispose.

Polyform-FTM may be used in conjunction with Formalex® (OPTIONAL EXTRA), a liquid control agent in a spray bottle, for small spills, vapour control for trash cans, and for hard to reach areas affected by the spill, such as cabinet faces and drawers. Polyform F neutralises 4-5 times its volume of formalin, formaldehyde and glutaraldehyde solutions.

1 - FAN PAD GLTM, Formaldehyde Neutralising Pad super absorbent pads and sheets, specially treated with Formalex® for the ultimate destruction of formaldehyde, 10% formalin and glutaraldehyde solutions and eliminate their hazardous vapors. Absorbs & neutralises up to 900cc of 10% formalin.



**Liquid Spill Lock
Formaldehyde Spill Kit**

Product Code: 258510FH

Qty	Items
1	Instructions
1	Compact carry bag
1	100ml Liquid Spill Lock formaldehyde neutraliser shaker
1	Fan pad GL
1	Pair protective gloves
1	Face mask with eye shield
1	Disposable apron
2	Disposable cloth wipes
2	Disposable waste bag & cable ties
2	Disposable scoop & scraper



Formaldehyde Neutralising Pads



Liquid Spill Lock Formaldehyde Neutraliser

Liquid Spill Lock - Mercury Spill Kit

A safer way to deal with Mercury spills

"Mercsorb" Amalgamation powder was developed to make cleaning mercury spills from work surfaces, floor cracks and other hard to reach areas easier and safer. This is accomplished through converting the elemental mercury into an amalgam, which ceases the emission of dangerous mercury vapors from being emitted.

This kit absorbs approx 25ml of Mercury.
Mercury Spill Kits are for single use only.

CAUTION - The contents of Mercury Spill Kits are not designed to make the mercury less toxic or dangerous.

Safety spill kit containing everything needed to contain, remove and decontaminate all stages of mercury spill cleanup. Provides a quick, effective and safe method of handling all sizes up to 25ml of mercury spills in labs, manufacturing facilities and anywhere else they may occur. Removes mercury from work surfaces including bench tops, floors and hard-to-reach areas such as cracks and crevices.

- Mineral, mercury absorbing powder
- "Mercsorb" Sponges also assist with clean-up
- Mercury Indicator powder assists with identifying any residual mercury
- Basic Personal Protection Equipment and clean up implements supplied

CAUTION - The use of protective breathing apparatus is recommended when dealing with Mercury Spills as the vapour is toxic. Please consult Silverback for supply of the appropriate item.



Liquid Spill Lock Mercury Spill Kit

Product Code: 258510MER

Qty	Items
1	Instructions
1	Compact carry box
1	500g Mercsorb amalgamation powder
1	Bottle water for powder mixing
1	150g Mercsorb indicator powder
2	Mercsorb sponges
1	Mercsorb waste container
1	Mixing tub with spatula
1	Pair protective gloves
1	Face mask with eye shield
1	Disposable apron
1	Disposable waste bag & cable tie
1	Disposable scoop & scraper



Mercsorb Sponges



Mercsorb Indicator Powder

LSL Compact Sharps Collection Kit

The name says it all.

The LSL Compact Sharps Collection Kit is designed for keeping you safe when removing unwanted sharps.

The fail safe picker obviates the need to handle syringes directly. The compact syringe disposal container is designed to capture the syringe with minimal risk of needle-stick injuries. The syringe container conforms to AS4931:1996 under manufacturing standards ISO 9001:2008.

A very inexpensive way of ensuring the safety of the user.



LSL Compact Sharps Collection Kit

Product Code: 25855SH

Qty	Items
1	"FITPACK" 3 Sharps disposable container
1	Sharps collection picker
1	Pair Nitrile disposable gloves
1	Absorbent cloth wipe
1	Clinical waste disposal bag with self seal strip
1	Instructions for use





Silverback Prencos Spill Kits & Services

Silverback offers an extensive range of Spill Kits suitable for use in a wide variety of environments such as medical kits featured in this brochure, Oil & Fuel, General Purpose and Hazardous Chemical kits.

Our Spill Kits can be offered in handy single use packs, compact and versatile portable kits presented in carry bags and buckets, larger wheelie bins for improved storage and mobility and towable trailers for large scale industrial applications.

Silverback also provides a range of services such as site inspections, Spill Kit customisation options, Spill Kit servicing and replenishment and Spill Kit/ Product training for employees and designated Safety Wardens.

For more information about our Spill Kit offering please contact your local Silverback Branch and speak to one of our friendly and helpful team members.



1300 858 858

www.prenco.com.au

www.silverback.com.au

A.B.N. 26 077 294 154

Melbourne

86 Strzelecki Avenue,
Sunshine West VIC 3020 Australia
T: 1300 858 858
E: melbourne@silverback.com.au

Sydney

26 Sleigh Place,
Wetherill Park NSW 2164 Australia
T: 1300 858 858
E: sydney@silverback.com.au

Brisbane

Unit 1, 12 Boolarra Street
Hemmant QLD 4174 Australia
T: 1300 858 858
E: brisbane@silverback.com.au

Adelaide

27 Pambula Street,
Regency Park SA 5010 Australia
T: 1300 858 858
E: adelaide@silverback.com.au

Auckland

300 Richmond Road,
Grey Lynn 1147 P.O.Box 46018
Auckland, New Zealand
T: +64 9 320 4090
E: auckland@silverback.co.nz

Head Office Admin

147 Railway Place,
Williamstown VIC 3016 Australia
T: +61 3 9312 0388